

An HMK-22 Fisher Sub-Sieve Sizer includes:

- HMK-22 host,
- Sample Pressure Calibrator,
- Silent air pump,
- Pressure regulator,
- Double-range flowmeter,
- Calculator chart,
- Precision-bore stainless steel sample tube,
- Two porous plugs,
- Plug manipulator,
- Plastic 2 oz. (60mL) funnel,
- Sample tube support,
- Drying agent,
- 3/4 in. diameter (19mm) paper discs and operating instructions.

[Back To Fisher Sub-Sieve Sizer](#)